# **QRT-1 SERIES**



MAX / DELUXE / STANDARD 4-POINT SECUREMENT WHEELCHAIR RETRACTORS



# **QRT-1 SERIES**

## The Securement System That Changed Everything

The original 4-point wheelchair securement system, QRT-1 Series retractors defined the way passenger safety devices are designed and tested.

With a range of 3 different options to fit every need and every budget, the QRT-1 Series of retractors offer easy to use, effective 4-Point securement of wheelchairs for virtually any vehicle application.

# QRT MAX

FULLY AUTOMATIC, premium knobless retractor that allows for one-handed operation.

## QRT DELUXE

SELF-LOCKING & SELF-TENSIONING retractor with ergonomic housing and dual tensioning knobs.



## QRT STANDARD

SEMI-AUTOMATIC retractor that meets all industry standards and specifications.

More than 30 years ago, Q'STRAINT introduced the world's first fully integrated 4-Point wheelchair passenger securement system, now an industry standard the world over. The QRT line of retractors are the linchpin of that system.





Every QRT retractor is fully ADA complaint, and meets or exceeds all standards and regulations, including:

- SAE J2249, ISO 10542,
- FMVSS 209, 302, 210, 222
- CMVSS 209
- CSA Z605
- and 30mph/20g crash testing

### Anchorage Options

All QRT-1 Series Retractors are compatible with L-Track, L-Pockets and Slide 'N Click anchorages, or may be directly mounted to vehicle floors, seat legs or barriers.







SLIDE 'N CLICK For kits that include Slide 'N Click anchorages, QRT Series retractors feature a single-bolt SNC assembly and plunger that allows a full 360° rotation, eliminating anchorage alignment guesswork.

## **QRT-1 SERIES** FEATURES COMPARISON

	ΜΑΧ	DLX	STD	
Knobless, One-Handed Operation. No knobs to interfere with wheels and footrests.	0			
Dual Tensioning Knobs. Provides additional tensioning if needed.		0		
Single Tensioning Knob. Provides additional tensioning if needed.			0	
Automatic, Self-Locking. Allows easy, one-handed hook-up.	0	0		
Self-Tensioning. Retractors automatically take up 'slack'.	0	0		
Positive Lock Indicator. Patented feature clearly indicates when fitting is locked in anchorage.	0	0	0	
Interchangeable. Eliminates confusion: no right, left, front or rear locations.	0	0	•	
Low Profile & Compact. Elimination of mounting bracket allows retractors to fit under most footrests.	0	0	0	
Accommodates Larger Wheelchairs. Reduced overall length leaves more room for wheelchairs.	0	0		
Ultra-Durable. Hardened steel and coated zinc for maximum corrosion resistance.	0	0	0	
Universal Design. Accommodates virtually all wheelchair designs, including scooters.	0	0	0	
J-Hook. Reduces twisting of belts and ensures proper securement for all wheelchair designs.	0	0	0	
Foot Release Lever. Easy release eliminates the stress of bending down.	0	0	0	



### WWW.QSTRAINT.COM/QRT-1-SERIES

© 2017 Q'STRAINT®. All Rights Reserved. Q'STRAINT®, Q'STRAINT® logo, QRT-SERIES® and QRT™ are registered trademarks of Q'STRAINT®, Inc.

#### **Q'STRAINT EUROPE**

72-76 John Wilson Business Park Whitstable, Kent, CT5 3QT United Kingdom Tel: +44 (0)1227 773035 Fax: +44 (0)1227 770035 Email: info@qstraint.co.uk

#### **Q'STRAINT AMERICA**

4031 NE 12th Terrace Oakland Park, FL 33334 Tel: 800-987-9987 Fax: 954-986-0021 Email: qstraint@qstraint.com

#### Q'STRAINT AUSTRALIA

Tramanco Pty Ltd. 21 Shoebury Street, Rocklea, Australia, QLD. 4106 Tel: +61 7 3892 2311 Fax: +61 7 3892 1819 Email: info@tramanco.com.au